

Quartz Screw-Top Sample Vials

MQSV-20-QCV | Features and requirements to define



1 Threaded neck and cap contact

Identify the mating neck, thread and closing contact.

2 Quartz bottle body

Specify the bottle envelope separately from the installed cap.

3 Cap removal space

Allow room to remove the cap and approach the sample.

REFERENCE VALUES

Body diameter

28 mm

Marked body height

58 mm

Illustrative design example; not an actual customer requirement.

Not to scale. Product-reference rendering with annotations; use the agreed drawing for dimensional acceptance.