


Micro Polarimeter Tube

MQOPT022 | Features and requirements to define



1 End-window interface

Confirm the end-window and sealing interfaces.

2 Small-bore sample region

Define the small-bore sample region.

3 Loading amount versus chamber value

Separate the required fill from the chamber value.

REFERENCE VALUES

Nominal path option

50 mm

Marked body diameters

OD14 / ID4 mm

Illustrative design example; not an actual customer requirement.

Not to scale. Product-reference rendering with annotations; use the agreed drawing for dimensional acceptance.